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Contents

1. Introduction	13
1.1. Foreword	13
1.2. Customer information	13
1.3. Mode of operation	13
1.4. Application	13
1.5. Scope of delivery	14
1.6. Declaration of conformity	14
1.7. Features	14
1.8. Indications	14
1.9. Contraindications	14
1.10. Safety Instructions	14
1.11. Warranty	15
2. Individual Adjustments	16
3. Cleaning, maintenance and disinfection	18
3.1. Care	18
4. Technical data	19
5. Choice of size	19
6. Size chart and article numbers	19
7. Transfer of the brace	20
8. Disposal	20

User Instructions

1. Introduction

1.1. Foreword

Our MKS[®] Modular Corset System straightens and stabilizes defined areas of the vertebral column in order to achieve postural corrections which actively support rehabilitation in specific indications.

Our MKS[®] programme is comprised of a modular system which can be divided into an orthosis set consisting of rigid supportive elements and a set of partially elastic, smoothly-fitting back-and-body bandages which adapt comfortably to the body.

By combining the individual components of the two sets an exact adjustment to both the respective indication and the patient can be achieved.

1.2. Customer information

For your own safety please read through these User Instructions carefully and accurately before using the brace. The instructions, notes and procedures must be read and understood thoroughly in order to benefit from the correct operation and use of the device. If anything in the User Instructions is not clear, or any instructions, operating procedures or safety information is not fully understandable, please contact the appropriate specialist retailer or albrecht GmbH directly, before you use the brace. This particularly applies to the safety instructions.

1.3. Mode of operation

The MKS[®] Pontsana[®] functions according to the bridging principle. The wide bandage with the additional lateral double straps straightens and hence relieves the lumbar vertebral column.

1.4. Application

The brace must only be used for the orthotic treatment of the vertebral column.

The MKS[®] Pontsana[®] assists treatment in the application area from the eleventh thoracic vertebral body to the fifth lumbar vertebral body. It provides stabilisation of the entire spine. The MKS[®] Pontsana[®] is intended to help you achieve a lifestyle that is as independent as possible, while still meeting your needs and matching your personal capabilities.

BWK 11 bis LWK 5

MKS® Pontsana

User Instructions

1.5. Scope of delivery

The MKS® Pontsana® is delivered to you ready to use in the size requested, with User Instructions and labeling on the product.

1.6. Declaration of conformity

The albrecht GmbH company, as the manufacturer solely responsible, declares that the MKS® Pontsana® brace conforms to the Regulation (EU) 2017/745 concerning medical devices.

1.7. Features

- Semi-elastic bandaging material
- Firm abdominal piece with hook and loop fastener
- Lateral reinforcing stays
- Bridging rods (spring band steel)
- Additional lateral double divided straps

1.8. Indications

The physician will decide on the type of treatment to apply based on his diagnosis. Generally, the use of the MKS® Pontsana® is indicated in all indications requiring a partial relief by delordozation of the lumbar spine, such as:

- Lower back pain
- Moderatley severe lumbar sciatica with slight failures in intervertebral disc protrusion/prolapse
- Moderatley severe radicular lumbar syndrome, pseudo-radicular lumbar syndrome

- Spondylolithesis, grade I with lumbago
- Moderatley severe lumbar deformity (relapsing treatment, complicated course)
- Condition after intervertebral disc surgery with small up to medium-sized prolapse
- Spondylarthrosis
- Sacroiliac joint irritation
- Lumbar muscle insufficiency

For all other indications a physician must be consulted.

1.9. Contraindications

The brace is only intended for use in contact with intact skin.

1.10. Safety Instructions

The optimal effect of the brace is only achieved when used correctly.

- The brace must only be used in the intact, complete and mechanically undamaged condition and with complete and intact cushioning and walers. This must be verified by the user before each usage.
- Opening or removing one or more belts, as well as excessive loosening of the waler when using the brace leads to a reduction of the therapeutic effect of the brace and may lead to injury.
- The brace must not be worn over open wounds.
- The skin should be free of oils, grease, gels or other debris, to prevent reactions with the skin or the structure of the material.

User Instructions

- The orthosis should fit firmly but not too tight, so as not to restrict the blood circulation and adversely affect nerve and lymph vessels. Excessive compression is therefore to be avoided.
- Combination with other products is currently not provided for or is to be agreed with the manufacturer in writing.
- The brace is not intended for single use, but is intended for multiple use by a single person.
- The product as delivered is not sterile.
- Contact your physician immediately in the event of an allergic reaction.
- Please note that cushioned sections can heat up under direct sunlight. Protect the orthosis from direct sunlight if necessary.
- Currently there is no test for flammability. Exercise caution when using the orthosis in the direct vicinity of open flames such as lighters and cigarettes.

1.11. Warranty

We assume the warranty for the brace for a period of 3 months. The brace is a medical rehabilitation device and must not be used for any purpose other than the intended, as described in the Instructions for Use. Changes to the brace or other applications require the written consent of the albrecht GmbH. If this is not obtained, the manufacturer may not honor the guarantee. If you are using single joints or other components, these should be used as intended. If changes or modifications (e.g. additional mounting holes) are made to the individual parts or components, the manufacturer's warranty no longer applies. The removal of or damage to the QM seal will also void the warranty.

2. Individual Adjustments

The MKS® Pontsana® is adjusted to the needs of patients and indications by adjusting the bridging rods by your orthopedic specialist.

For those patients using the MKS® Pontsana® for the first time it is important that they learn how to handle it. Putting it on and taking it off should be practiced.

Wear the MKS® Pontsana® directly on the skin or over a tight-fitting singlet without any creases. It is recommended that an absorbent singlet be worn to avoid humidity caused by excessive sweating.

To put on the MKS® Pontsana® please follow the instructions below. This will ensure reliable effectiveness.

Always make sure that the MKS® Pontsana® is a smooth and tight fit. However, in order to prevent constrictions along the edges it shouldn't be over-tightened.

User Instructions

Open the inner and the outer hook and loop fastener. Put the bandage around your waist so that the openings of the pockets for the stabilising rods point upwards.

- 1 Press the left-hand side fluffy band of the inner bandage belt flat on your stomach, and then pull the hook and loop end of the right-hand side inner strap forward and attach it.
- 2 Pull the outer strap ends forward to increase the stabilising pressure around your waist.
- 3 Press the tightened outer strap ends flat onto the closed inner strap and flatten it down with your flat hand.

Hand wash at 30°C

Do not bleach

Do not iron

Do not dry-clean

Do not tumble dry

3. Cleaning, maintenance and disinfection

The orthosis is designed to be maintenance-free. To ensure proper operation over the period of treatment the orthosis should be cleaned regularly (at least every 3 months) or as required, according to the following instructions.

3.1. Care

- The product as delivered is not sterile and is not suitable for sterilization.
- Do not use harsh or abrasive cleaners.
- All fabrics can be washed by hand at 30°C using water and a mild detergent and/or disinfectant.
- Not machine washable.
- The MKS® Pontsana® can be dried in the open air.
- In the case of more severe soiling, a replacement set of textile parts is available.
- Wipe down surfaces with a cloth soaked with disinfectant.
- Wet completely, and do not wipe off.
- Spray inaccessible surfaces.
- When spraying ensure complete wetting.
- A mild alcohol-based disinfectant is recommended.

Ask your physician or pharmacist when selecting a disinfectant, and follow the instructions given by the disinfectant manufacturer. The Robert Koch list of approved disinfectants can be found at <http://www.rki.de>.

User Instructions

4. Technical data

Description	Material
	elastic material 28% viscose, 48% perlon, 24% elastodiene
Body Belt	binding 83% viscose, 17% dorlastan
	velours PA

5. Choice of size

The MKS[®] Pontsana[®] is available in sizes according to the table below.

6. Size chart and article numbers

Size	Body circumference	Abdominal height	Height of back piece	Art.-No.
1	60 - 70 cm	18 cm	32 cm	582-801-S
2	70 - 80 cm	18 cm	32 cm	582-802-S
3	80 - 90 cm	18 cm	34 cm	582-803-S
4	90 - 105 cm	18 cm	34 cm	582-804-S
5	105 - 120 cm	18 cm	38 cm	582-805-S
6	120 - 135 cm	18 cm	38 cm	582-806-S
51	105 - 120 cm	18 cm	34 cm (reduced)	582-8051-S
61	120 - 135 cm	18 cm	34 cm (reduced)	582-8061-S

MKS[®] Pontsana[®] – reduced height:

Size	Body circumference	Abdominal height	Height of back piece	Art.-Nr.:
1	60 - 70 cm	16 cm	26 cm	582-801N-S
2	70 - 80 cm	16 cm	26 cm	582-802N-S
3	80 - 90 cm	16 cm	26 cm	582-803N-S
4	90 - 105 cm	16 cm	26 cm	582-804N-S
5	105 - 120 cm	16 cm	26 cm	582-805N-S
6	120 - 135 cm	16 cm	26 cm	582-806N-S

7. Transfer of the brace

The brace is not intended for single use, but rather is intended for multiple use by a single person. We do not recommend transfer to other users. Should this be desired however, please ensure to pass on the care and cleaning instructions and have the hinge checked by an authorized specialist dealer for safe and proper operation.

8. Disposal

The brace contains recyclable materials without toxic or other harmful substances or other environmentally hazardous substances. Provided it is not contaminated with infectious germs, the hinge can be deposited in the normal waste disposal. To be sure, consult your specialist orthopedics dealer.

Duty to report

Due to regional legal regulations, you are required to immediately report any serious incident involving the use of this medical device to the manufacturer and the responsible authorities. Please find our contact details on the back of this brochure.

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